Botany Evolution LLC

CERTIFICATE OF ANALYSIS

Report Date	9-Aug	Country of Origin	Sol	Solomon Islands	
Sample Number	S2164	Country of Processing		USA	
Product Name Gol	d Basal Roots	Manufacture Date		Jun-24	
Lot Number SI	<2406-GBR8	Best By Date		Jun-27	
ITEM	SPECIFICATION	TEST RESULT	S	METHOD	
IYSICAL & CHEMICAL					
Identification	Piper methysticum	Complies		HPLC	
Appearance	Beige to Yellow Powder	Complies		Organoleptic	
Kavalactone Standard	2-17 % Kavalactones	9.90%	And a	HPLC	
Kavalactone Profile	Noble	Pass	All and	HPLC	
Chemotype	If # 5 is in 1st or 2nd in Abundance	243516	A AND	HPLC	
K/DHM	> 1.2 for Noble	4.0	9.69/14	Calculation	
EAVY METALS				5 77	
		Results		CP120	
Arsenic (As)	NMT 1,000 (ppb)*	< 10	ppb	FDA EAM 4.7	
Cadmium (Cd)	NMT 1,000 (ppb)*	190	ppb	FDA EAM 4.7	
Lead (Pb)	NMT 1,000 (ppb)*	1,050	ppb	FDA EAM 4.7	
Mercury (Hg)	NMT 1,000 (ppb)*	< 10 ppb		FDA EAM 4.7	
*Heavy Metals Action Limits Based on I	Maximum PDE at 5% Kavalactones. Results M	ay Exceed 1,000 ppb action li	mit with higher l	kavalactone contents.	
ICROBIOLOGICAL	1 67				
ICROBIOLOGICAL		Results			
AEROBIC PLATE COUNT	NMT 10,000,000 cfu	11,000	cfu / 10 g	USP 2021	
E. COLI	ABSENT (cfu/10g)	Absent	cfu / 10 g	USP 2022	
PSEUDOMONAS AERUGINOSA	ABSENT (cfu/10g)	Absent	cfu / 10 g	USP 2022	
SALMONELLA	ABSENT (cfu/10g)	Absent	cfu/10g	USP 2022	
STAPHYLOCOCCUS AUREUS	ABSENT (cfu/10g)	Absent	cfu / 10 g	USP 2022	
YEAST		8,000	cfu / 10 g	0	
MOLD	NMT 100,000 cfu (Combined)	240	cfu / 10 g	USP 2021	
TOTAL YEAST & MOLD	NMT 100,000 cfu (Combined)	8240	cfu / 10 g		
		al de recent			

Analysis Performed by a Third-Party Laboratory

We are dedicated to offer the best quality of botanical products on the market. We test and stand behind our products. Disclaimer * The test results are accurate to the best of our knowledge and are based upon reputable laboratory and industry standard testing methods.

These results should not be used as a final determination for use in a finished product. It is recommended that you verify these test results with an in-house quality control department or obtain an additional independent third-party lab to verify that this material meets specifications.

Botany Evolution, its board of directors, contract laboratories, employees, and affiliates are held harmless from any loss or damages resulting from the use or misuse of this document. The appropriate use of this product is the sole responsibility of the user of the purchasing party.

Completed By: ////

Botany Evolution LLC

2510 Kirby Circle NE alm Bay, FL 32945 21-802-4583

Certificate Of Analysis

Sample Identification Information

		and the second second second second						
Date of Analysis 8/9/2024				Country of	Origin	SOLOMON ISLANDS		
Sample: S2164			<u>(</u>	Country of Pro	ocessing	USA		
Product Name GOLD BASAL				Manufacture Date		Jun-24		
Lot# SIK2406-GBR8				Best By Date		Jun-27		
General Product Spec	ificatio	ons						
Product Species Piper Methysticum				Common Names		Kava kava, Awa, Yagona		
Part Used Root				Appearance		Yellow, Brown, beige powder		
Analyzed Characteris	tics	Specificat	ion		<u>Result</u>	1 min	Test Method	
Standardization		2-17% Kaval	actones		9.90%	" (150)	HPLC	
Identification		Complies by	HPLC		Conform		HPLC	
Kavalactone Profile		Noble			PASS		HPLC	
Mesh Size	(\bigcirc)	60-30		D/	60		Sieve	
<u>Color</u>		Beige to Yell	ow		Pass		Visual	
<u>Odor</u>			1.6		Pass		Organoleptic	
<u>Taste</u>				11	Pass		Organoleptic	
<u>Chemotype</u>					243516		HPLC	
<u>K/DHM</u>			1		4.0		Calculation	
Kavalactones	Code	Peaks Ref.	Correction	Area *	Area %	Corrected	Chemotype Identi	

Kavalactones	Code	Peaks Ref. (elution order)	Correction Factor	Area *	Area %	Corrected Kavalactones	Chemotype Identifier
Standard Kavain	к			2372.192	1.0		
Methysticin	M	1	2.21	350.765	3.93%	0.45%	6
Dihydromethysticin	DHM	2	3.38	389.123	4.36%	0.77%	5
Kavain	K	3	1	5221.783	58.57%	3.06%	4
Dihydrokavain	DHK	4	3.48	1711.079	19.19%	3.49%	2
Desmethoxyyangonin	DMY	5	2.52	448.729	5.03%	0.66%	1
Yangonin	Y	6	3.12	793.475	8.90%	1.45%	3
Kavalactones	6		Total:	8914.954	100.00%	9.90%	243516

8/13/24

Date

*See data in attachment HPLC1100 Agilent Certificate with Chromatogram graph.

This result are in house tested and the best of our knowledge and experience. Using calibrated equipment.

Mustl

We are dedicated to offer the best Quality of Botanical products on the market. We test and stand behind our products.

Disclaimer* the test results are accurate to the best of our knowledge and are based upon reputable laboratory and industry standard testing methods.

These results should not be used as a final determination for use in a finished product. It is recommended that you verify these test results with an

in house quality control department or obtain an additional independent third party lab to verify that this material meets specifications

Botany Evolution, its board of directors, contract laboratories, employees, and affiliates are held harmless from any loss or damages resulting from the

use or misuse of this document. The appropriate use of this product is the sole responsibility of the user of the purchasing party.

Yleings

Chemist

otany Evolution LLC	
510 Kirby Circle NE	Kavalactone Analysis
alm Bay, FL 32945	
21-802-4583	SAMPLE S2164
	Vial 11
.75325g/50mL	
avelength 246 nm	
	08_09_2024_15MINSTDTESTMETHOD 2024-08-09 12-47-50\01->
EQUENCE C:\CHEM32\1\DA	
njection date 8/9/202	
njection time 4:34:47	PM
.cq. operator KRISTL	
lethod	C:\CHEM32\1\DATA\KAVA_08_09_202->
	=500,60 (KAVA_08_09_2024_15MINSTDTESTMETHOD 2024-08-09 12-47-50\011-1101.D)
mAU	6.404 - kavain
400	ž.
350	8
	6.40
300	DAR C
250	E S S
250	dihydrokavain
200	
	- methysticin - dihydromethysticin - desmethoxyangonir
150	and an an and an
100 -	
	9.194
50	
0	
0 2.5	5 7.5 10 12.5 min
# COMPOUND	RET. TIME AREA AREA & AMOUNT
1 methysticin	5.562 350.765 3.93 0.000
2 dihydromethysticin	
3 kavain 4 dihydrokavain	6.4045221.78358.570.0007.0131711.07919.190.000
5 desmethoxyyangonin	8.925 448.729 5.03 0.000
6 yangonin	9.194 793.475 8.90 0.000
5 yangonin	9.194 193.475 8.90 0.000
	S S A A A A A A A A A A A A A A A A A A
말 승규는 것 같아요. 이는 것 같아.	

8/13/24